

MICROPROCESSOR CONDUCTIVITY / TDS / TEMP.

METZ – 1001 M



SPECIFICATIONS

Display	16 X 2 Alpha Numeric LCD Display.
Range	0.1 μ S to 100 ms (6 Decadic Ranges).
Resolution	0.01 μ S.
Accuracy	$\pm 1\%$ of F.S. ± 1 Digit.
Calibration	Auto / Manual.
Temperature Range	0 to 100°C.
Cell Constant	0.1 to 5.0 Adjustable.
Storage	Up to 100 Samples.
Power	230V $\pm$ 50Hz AC.

METZER OPTICAL INSTRUMENTS PVT. LTD.

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

www.metzterlab.com info@metzterlab.com +91 9897000955