

MICROPROCESSOR BASED pH METER

METZ – 2001 M



SPECIFICATIONS

Display	16 X 2 Alpha Numeric LCD Display
Range	pH: 0 to 14.00, mV : 0 to ± 1999.9 .
Resolution	0.01 pH, 0.1 mV, Temp. 0.1°C
Accuracy	± 0.01 pH, $\pm .1$ mV, Temp. $\pm 1^\circ\text{C}$. (Auto/Manual)
Temp. Compensation	0 to 100°C.
Calibration	Auto / Manual.
Storage	Up to 100 Samples.
Power Supply	230V $\pm 10\%$, 50Hz AC.
Auto Buffer Recognition	Auto 4.00, 7.00, 9.20 pH.

METZER OPTICAL INSTRUMENTS PVT. LTD.

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

🌐 www.metzlab.com 📧 info@metzerlab.com ☎ +91 9897000955