

AUTO DELUXE pH METER

METZ – 201 M



SPECIFICATIONS

Display	3.5 Digit LED.
Range	pH 0 to 14.00, mV 0 to ± 1999 . Touch Button Calibration.
Resolution	pH 0.01, mV 1.
Accuracy	pH ± 0.01 , mV ± 1 mV.
Temp. Compensation	Auto 0 to 100 deg. C., Manual 0 to 100 deg. C.
Slope Control	80 to 120%.
Power Supply	230V $\pm 10\%$ 50 Hz AC.

METZER OPTICAL INSTRUMENTS PVT. LTD.

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

🌐 www.metzerlab.com 📧 info@metzerlab.com ☎ +91 9897000955