

DIGITAL SPECTROPHOTOMETER (DOUBLE DISPLAY)

METZ – 3001 M



SPECIFICATIONS

Range	340 – 960nm.
Resolution	0.1 % T, 0.001 Abs.
Output	%T: 0 to 100%, Abs: 0 to 1.999.
Detector	Silicon Photodiode.
Light Source	6V – 1Amp. Halogen Lamp.
Wavelength Resolution	1 nm.
Repeatability	± 1nm.
Display	3.5 Digit LED Display.
Power Supply	230V ± 10% 50 Hz AC.

METZER OPTICAL INSTRUMENTS PVT. LTD.

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

www.metzlab.com info@metzlab.com +91 9897000955