

SINGLE BEAM UV-VIS SPECTROPHOTOMETER

METZ – 3003 M



SPECIFICATIONS

Photometric Range	% T:0 TO 100%, Abs:- 0 to 1.999, Conce, 0 to 1.999 K-Factor: 0 to 1.999.
Accuracy	±0.005 Abs at 1.0 Abs ±0.010 Abs at 1.0 Abs.
Data Readout	Single Line LCD display.
Key Board	4 Keys, Soft touch membrane type.
Light Source	Tungsten Halogen Lamp.
Detector	Wide range Silicon Photodiode.
Sample Holder	4 position adjustable sample holder.
Wavelength	Selenium Photocell / Photodiode.
Light Source	195 to 1000nm.
Spectral Bandwidth	2 nm.

METZER OPTICAL INSTRUMENTS PVT. LTD.

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

www.metzlab.com info@metzlab.com +91 9897000955