

DOUBLE BEAM UV-VIS SPECTROPHOTOMETER

METZ – 3004 M



SPECIFICATIONS

Optical System	Double Beam , Grating 1200 lines / MM.
Wavelength Range	190-1100nm.
Bandwidth & Wavelength Accuracy	1nm & 0.3 nm.
Wavelength Scanning	Auto.
Key Board	Membrane keypad.
Light Source	Deuterium & Tungsten lamp.
Display	320 x 240 LCD with graphic Display.
Detector	Silicon Photodiode.
Power supply	AC 220V/50 Hz.

METZER OPTICAL INSTRUMENTS PVT. LTD.

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

🌐 www.metzlab.com 📧 info@metzerlab.com ☎ +91 9897000955