

DIGITAL PHOTO COLORIMETER WITH AUTO ZERO CALIBRATION

METZ – 401 M



SPECIFICATIONS

Standard Glass Filters	8 (Wavelength range 400 – 700nm).
Minimum Volume	1ml.
Display	3.5Digit Bright Red 7 Segment LED Display
Key Board	1 Keys Soft Touch Membrane Type.
Absorbance	0 – 1.99.
Resolution	0.01
Accuracy	0.5 % FSI.
Detector	Selenium Photocell / Photodiode.
Light Source	6.3V – 0.3 Amp. Tungsten Lamp.
Power Supply	230V ± 10%, 50Hz. AC.

METZER OPTICAL INSTRUMENTS PVT. LTD.

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

🌐 www.metzlab.com 📧 info@metzerlab.com ☎ +91 9897000955