

M DIGITAL FLAME PHOTO METER (DUAL CHANNEL Na & K)

METZ – 802 M



SPECIFICATIONS

Range	Na : 0 – 100 ppm & K : 0 –100 ppm.
Sensitivity	Na : 5 ppm & K : 5 ppm
Display	Dual 2 ½ digit bright LED Display
Accuracy	± 2% upto 100 meq/l , ± 5% above 100 meq/l
Repeatability	± 2 Count
Detector	Silicon Photodiode.
Filters	Narrow band interference glass filters.
Flame System	LPG and dry oil free air.
Warm up time	10 minutes.
Power Supply	230V ± 10% AC 50 Hz

METZER OPTICAL INSTRUMENTS PVT. LTD.

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

www.metzerlab.com
info@metzerlab.com
[+91 9897000955](tel:+919897000955)