

MICRO PROCESSOR BASED FLAME PHOTO METER

METZ – 803 M



SPECIFICATIONS

Range & Resolution	Na : 0.1 – 100 ppm & K : 0.1 –100 ppm& 0.1 ppm/meq
Sensitivity	Na : 0.5 ppm & K : 0.5 ppm
Display	Graphical LCD Display
Accuracy	± 2% upto 100 meq/l , ± 5% above 100 meq/l
Reproducibility	2% for 25 sample
Detector	Silicon Photodiode.
Filters	Filters for Sodium and Potassium.
Flame System	LPG and dry oil free air.
Range Selection	Automatic
Calibration	Up to 5-point calibration with curve fitting software
Power Supply	230V ± 10% AC, 50/60 Hz

METZER OPTICAL INSTRUMENTS PVT. LTD.

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

www.metzlab.com info@metzerlab.com +91 9897000955