

DIGITAL DISSOLVED OXYGEN METER

METZ – 902 M



SPECIFICATIONS

Display	3.5 Digit LED.
Range	0 to 20.0ppm
Resolution	0.1ppm
Accuracy	± 0.2ppm
Temp. Compensation	0 to 50 deg.C.
Sensor	DO Amperometric Gold / Silver
Time Response	10 sec.
Power Supply	230V ± 10% 50 Hz AC

METZER OPTICAL INSTRUMENTS PVT. LTD.

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

🌐 www.metzerlab.com 📧 info@metzerlab.com ☎ +91 9897000955