

DIGITAL HAEMOGLOBIN METER

METZ – 903 M



SPECIFICATIONS

Display	Digital 3½ Digit, Bright LED
Display Mode	% T , O.D. and 3 Concentration Modes C1 , C2 , C3
Zero Adjustment	Automatic (Without Knob Rotation)
Accuracy	± 0.02 O.D.
Light Source	6.8 Volt, 0.3 Amp Tungsten Lamp
Memory	In built memory to retain Last "K" factor factor for concentration even after power loss.
Sample Volume	1 ml. Minimum.
Power Supply	230V ± 10% AC 50 Hz

METZER OPTICAL INSTRUMENTS PVT. LTD.

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

www.metzlab.com info@metzlab.com +91 9897000955