

DIGITAL DRY BATH INCUBATOR

METZ – 904 M



SPECIFICATIONS

Display	LCD Display for time & temperature.
Block size	24 Holes for 12mm tubes
Temperature Range	+5 °C above ambient to 100°C
Temperature Accuracy	± 2°C
Maximum heating power	125 W
Timer	1 min to 99 min.
Power Supply	230 V ± 10% AC , 50Hz

METZER OPTICAL INSTRUMENTS PVT. LTD.

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

🌐 www.metzlab.com 📧 info@metzerlab.com ☎ +91 9897000955