

# AUTO DIGITAL pH, CONDUCTIVITY & TEMP. METER

## METZ –203 M



### SPECIFICATIONS

Display	3.5 Digit LED.
Power	230V ± 10% 50 Hz AC.
PH :	
Range & Resolution	pH 0 to 14.00 & pH 0.01.
Accuracy & Temp. Comp.	pH 0.01 ± 1 Digit & 0 to 100 Deg. C.
Conductivity :	
Range	0 to 20mS / cm in 2 ranges.
Resolution	0.001 m Mhos (1Micro Mhos).
Accuracy & Cell Constant	0.5% of range ± 1 digit & adjustable.
Temperature :	
Range & Resolution	0 to 100 deg.C & 0.1 Deg. C.
Accuracy	± 0.2% of range ± 1 digit.
Sensor	PTD (PT-100) Probe.

**METZER OPTICAL INSTRUMENTS PVT. LTD.**

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

[www.metzerlab.com](http://www.metzerlab.com)
[info@metzerlab.com](mailto:info@metzerlab.com)
[+91 9897000955](tel:+919897000955)