

DIGITAL TURBIDITY METER

METZ –501 M



SPECIFICATIONS

Display	3.5 Digit LED.
Power Supply	230V $\pm$ 10% 50 Hz.
Range	0 to 1000 NTU in 2 ranges.
Resolution	1 NTU.
Accuracy	$\pm$ 3% FS $\pm$ 1 Digit.
Sample System	30mm Clear Glass Test Tubes.
Light Source	6V – 1Amp. Tungsten Lamp.
Detector	Photocell / Photodiode.

METZER OPTICAL INSTRUMENTS PVT. LTD.

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

🌐 www.metzlab.com 📧 info@metzerlab.com ☎ +91 9897000955