

DIGITAL IMAGING CAMERA

METZ – 4K ULTRA HD



Features

- METZ-4K ULTRA HD is the ultimate live do play microscope camera for bright field image. DIGI 4K offer full, ultra-high definition 3840 X 2160 output (up to 30fps) with 1400 TV lines of resolution. This fast frame rate enables the operator to move samples without smear lag or filter. It automatically adjusts to changes in light intensity.
- 3840 x 2160 – Highest Resolution HD images. (1400 TV lines of resolution).
- Easy to use-1 button white balance. Factory preprogrammed menu settings for optimal results.

SPECIFICATION

Sensor size	Sensor 1/1.9" SONY
Pixels Size	1200µm
Frame rate	60 FPS
Resolution	4032 x 3024P
Video Resolution	1920 x 1080P
Gain Control	Support
Pixel Depth	12bit
Definition	4K UHD

Exposure	Manual/AUTO
White Balance	Manual/AUTO
Photo & Video	Support
Image format	TIF, BMP, JPG, RAW
Operating system	Windows 7/8/10/11 32/64bit
Camera adapter	C-mount
Power supply	12V 2A

METZER OPTICAL INSTRUMENTS PVT. LTD.

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

🌐 www.metzlab.com 📧 info@metzlab.com ☎ +91 9897000955