

DIGITAL IMAGING CAMERA

METZ – M DIGITAL TABLET



Features

- M- DIGITAL TABLET is the smart digital microscope tablet camera equipped with 5MP color CMOS Camera and 9.7 Inches Tablet.
- It is developed specifically for microscopy applications.

SPECIFICATION

Sensor	5MP color CMOS Camera sensor, Size 1/2.5 inch (4:3)	Video and Frame rate	640x480 @ 30fps, 1280x720@15fps
Effective Resolution	2592Horizontal x 1944Vertical	Software Function	Image, Video
Pixels Size	2.2µm × 2.2µm	Optional Output	HDMI
Color Depth	12-bit, on-chip	Power supply	DC 5V 2A
SD Card	Support	PAD Size	220*170*20mm
Shutter Type	Electronic	Tablet weight	1 kg
Exposure	Auto Exposure		

METZER OPTICAL INSTRUMENTS PVT. LTD.

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

www.metzlab.com info@metzerlab.com [+91 9897000955](tel:+919897000955)