

DIGITAL IMAGING CAMERA

METZ – M 5.0 PRO



SPECIFICATION

| | | |
|-------------------|---|--|
| Image Sensor | 5.1 Megapixel 2/3 Interline Sony ICX252 | |
| Frame Rate | 30fps at 640x480 | Features <ul style="list-style-type: none"> • Interfaced to a computer via USB 2.0 technology for rapid image download. • Plug-to-use, suitable with any standard biological, stereo and metallurgical microscope. • Ease of operation; with any Metzer – M Digital Imaging software • Standard C- Mount makes connection to any microscope easy. |
| Resolution | 2592 x 1944 | |
| Color Depth | 8 and 16-bit | |
| Power Requirement | USB | |
| Auto Exposure | Automatic Manual | |
| White Balance | Automatic Manual | |
| PC Interface | Interface Connector Standard USB 2.0 high - speed interface | |