


INVERTED METALLURGICAL MICROSCOPE

METZ – 5000 ITM




SPECIFICATIONS

| | |
|------------------|--|
| Viewing Head | Trinocular Observation head conveniently inclined at 45°. |
| Magnification | 50x – 450x |
| Mechanical Stage | Large graduated mechanical stage with low drive co – axial controls for X – Y movements on ball slides. |
| Focusing | Co – axial coarse & fine focusing module with slides travelling on ball bearing guide ways. |
| Illumination | Incident light through Epi – illuminator 12V – 50W Halogen lamp with field aperture diaphragm and filter slot. Continuously variable luminosity control through built in electronic transformer. |
| Objectives | Achromatic flat field metallurgical corrected and anti – mount fungus M5x, M10x, M20x & M45x. |
| Eye pieces | Compensating High eye point wide field eye piece 10x. |
| Filters | Green & Blue. |

OPTIONAL ACCESSORIES



DIGITAL IMAGING CAMERA



METZER – M INDUSTRY
PRO SOFTWARE

METZER OPTICAL INSTRUMENTS PVT. LTD.

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

 www.metzlab.com
 info@metzerlab.com
 +91 9897000955