



# MEDICAL RESEARCH MICROSCOPE

## METZ-38



### SPECIFICATIONS

Magnification	100x – 1000x.
Stand	Robust and stable inclinable through 90°, Mechanical tube length 160mm.
Focusing	Separate knobs for coarse and fine motion.
Stage	Fixed square stage size 120 x 120 mm. with co-axial graduated Mechanical stage for X & Y movements up to 55mm to 75mm.
Illumination	Plano Concave mirror in fork mount replaceable by Sub stage microscope lamp.
Condenser	Sub stage Abbe filter holder N.A. 1.25 condenser with Iris diaphragm and swing out filter holder focusable by rack & pinion.
Nose Piece	Dust proof triple revolving nose piece with positive click stops.
Objectives	Achromatic 10x, 45x & 100x oil (SL.).
Eye Pieces	W.F. 10x Anti Fungus Coated.

**METZER OPTICAL INSTRUMENTS PVT. LTD.**

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

 [www.metzlab.com](http://www.metzlab.com)
 [info@metzerlab.com](mailto:info@metzerlab.com)
 +91 9897000955