


PARTICLE ANALYSIS SOFTWARE



SPECIFICATIONS

Live Measurements:

Distance, Length, Width, Perimeter, Angle, Three Point Radius.

Count & Classification: Identification of Objects in an image, count them, obtain Several features Measurements. Objects identification by user or automatically. User Defined classification on basis of Size or intensity.

Particle Size Analysis: Various calculation & measurements available for selected Particle are as Follows: Dimensions, Area, Perimeter, Ferrite Length, Min/Max Radius, Thread Length, Thread Width, Fiber Length, Fiber Width.

Densitometry Analysis: Manual, Auto bright and auto dark methods to identity Intensity range Defined object to be measured. Roundness, Shape Orientation, Elongation, Equal Circular Di - Ameter, Equal Sphere Volume.

Report Generation: Direct printout with original image processed image & Tabular Results or Export To MS Office or Excel for further modification.

METZER OPTO INC.

Tower A-0007, ATS Bouquet, Sector - 132, NOIDA - 201305 (U.P)

www.metzzerlab.com sales@metzermicroscope.com +91 9897000955



Scan to download
complete catalogue