

METZER MICRO-MEASURE



SPECIFICATIONS

- **MEASUREMENT:** Spatial calibration, Line Measurements for Distance, Length, Width, Perimeter, Angle. Three point radius. Area by enclosed line controlled by four arrow keys available on keyboard arrows with zoomed preview. Measurements can be performed on live images.
- **REPORT GENERATION:** Direct printout with original image processed image & tabular results or export to MS office or Excel for further modification. Calibration scale included.

METZER OPTO INC.

Tower A-0007, ATS Bouquet, Sector - 132, NOIDA - 201305 (U.P)

www.metzzerlab.com sales@metzermicroscope.com +91 9897000955



Scan to download
complete catalogue