

TRINOCULAR STEREO ZOOM MICROSCOPE

METZ – 7000 TZM WITH BOOM STAND



SPECIFICATION

Magnification Range	7 to 45x (Standard) can be increased to 180x depending on different combination of eye piece & Auxiliary lens used.
Observation Head	Trinocular, 1x magnification with interpupillary and diopter adjustment.
Stand	Double-arm heavy duty boom for stereo microscopes allowing you to point microscope head in ANY direction desired.
Base	Heavy Cast steel base.
Focusing System	Rack & pinion arrangement coarse focusing with pillar clamp. Coarse adjustment app. 50mm.
Illumination	LED ring light with intensity controller.
Objectives Range	Zoom 0.7 to 4.5x.
Eye Pieces	10x paired. Extra Wide Field, Anti-fungal & Anti reflected.
Weight	20 kgs Approx.