

# TRINOCULAR STEREO ZOOM MICROSCOPE

METZ – 7000 TZM



## SPECIFICATION

Magnification Range	7 to 45x (Standard) can be increased to 180x depending on different combination of eye piece & Auxiliary used.
Observation Head	Trinocular, 1x magnification with interpupillary and diopter adjustment.
Stand	Pillar plain stand, sturdy & durable.
Focusing	Rack & pinion arrangement coarse focusing with pillar clamp. Coarse adjustment app. 50mm.
Illumination	Built in base LED 3W – 3V transmitted & reflected light with light intensity control system.
Objectives	Zoom 0.7 to 4.5x.
Eye Pieces	10x paired. Extra Wide Field, Anti –Fungal & Anti reflected.

**METZER OPTICAL INSTRUMENTS PVT. LTD.**

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

🌐 [www.metzlab.com](http://www.metzlab.com)    📧 [info@metzlab.com](mailto:info@metzlab.com)    ☎ +91 9897000955