

DISSECTING STEREO BINOCULAR MICROSCOPE

METZ – 211



SPECIFICATIONS

Magnification	20x – 40x
Observation	Binocular inclined at 45°, Rotatable 360°.
Stand	Flat duly balanced with milky glass plate and clips for specimens.
Focusing	Rack & Pinion with extra – large knobs.
Objectives	2x and 4x mounted on inter changeable slide.
Eye Pieces	Wide field 10x

METZER OPTICAL INSTRUMENTS PVT. LTD.

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

🌐 www.metzlab.com 📧 info@metzerlab.com ☎ +91 9897000955