

DISSECTING STEREO BINOCULAR MICROSCOPE

METZ – 214



SPECIFICATION:

Magnification	20x – 40x
Observation	Binocular inclined at 45°, Rotatable 360°.
Stand	Heavy and stable construction.
Focusing	Rack & Pinion with extra – large Knobs and slip clutch Mechanism with safety stop.
Illumination	Stand provides both incident and transmitted light with all electrical and transformers built inside base having 6V – 20W lamps.
Objectives	2x and 4x paired mounted on revolving turret, changing from low to high by simple rotation and no exchange slides.
Eye Pieces	Wide field 10x

METZER OPTICAL INSTRUMENTS PVT. LTD.

C-21, GURU NANAK NAGAR, JUNCTION ROAD, MATHURA – 281001

www.metzlab.com info@metzerlab.com [+91 9897000955](tel:+919897000955)